



Micron offers *Cxl* series fulfills all features of a Modern microscopes with latest techniques, Economical Design, fulfills the complete features of a Biological Microscopes with Three models i.e. MONO *Cxl*, BINO *Cxl* & TRINO *Cxl*.

Technical Specifications		
Viewing Tubes	BINO <i>Cxl</i>	Binocular head 45° inclined
	MONO <i>Cxl</i>	Monocular head 45° inclined
	TRINO <i>Cxl</i>	Monocular head 45° inclined
Nose Piece	Quadruple Ball bearing Nose piece.	
Mechanical Stage	Low Drive Co-axial control, having X & Y movement of 55 mm & 75 mm.	
Focussing System	Co-axial Coarse & Fine focussing with focussing lock.	
Condensor	Abbe Condensor N.A. 1.25, with iris diaphragm movable on Rack & Pinion.	
Illumination	LED 3watt with intensity control, input 220V. Battery backup (optional) or 6V-20W, halogen bulb (on demand)	
Objective (Ach)	4x, 10x, 40x S/L & 100x Oil (P.D.:37MM)	
Eye Piece	Wide Field 10x	
Accessories	Dust Cover, Halogen Lamp	
Packing	Duly Packed in Styrofoam Box	



MONO *Cxl*

The manufacturer reserves the right to make changes in instrument design in accordance with scientific and mechanical progress without notice and without obligation. All Micron Optik Instruments are warranted against defect in materials and workmanship for one year .

micron[®]
OPTIK
MICROSCOPES.